**IVD DNlite uPTM-FetA ELISA Kit**

REF: 8103106

**Enzyme immunoassay for the quantitative determination of human unique Fetuin-A with specific post translational modification (PTM) in urine**
**INTENDED USE**

The DNlite uPTM-FetA ELISA Kit is a colorimetric immunoassay to quantify unique Fetuin-A with specific post translational modification (hereafter called as E103) in urine. This product is a non-automated IVD for prognosis use and should be performed at qualified clinical laboratories by certified medical professionals, such as Medical Technologists. The concentration of E103 should be corrected by urine creatinine before applying to the clinic. The DNlite uPTM-FetA ELISA Kit is to be used in conjunction with clinical evaluation as an aid in assessing the prognosis of kidney function in type 2 diabetes and kidney transplantation recipients.

**INTENDED PATIENT POPULATION**

For patients with type 2 diabetes:	For kidney transplantation recipients:
<ul style="list-style-type: none"> <li>• Diagnosis of type 2 diabetes in accordance with American Diabetes Association</li> <li>• Age <math>\geq 18</math> years</li> <li>• Not dependency of renal replacement therapy</li> <li>• Without current known malignancy</li> <li>• Without current acute or chronic infections</li> <li>• Without primary or secondary renal diseases (other than DKD)</li> <li>• Without illicit drug use</li> </ul>	<ul style="list-style-type: none"> <li>• Age <math>\geq 18</math> years kidney transplant recipient with functioning graft beyond 1 year</li> <li>• Not dependency of renal replacement therapy</li> <li>• Without current known malignancy</li> <li>• Without current acute or chronic infections</li> </ul>

**PRINCIPLE OF THE TEST**

The DNlite uPTM-FetA ELISA Kit is a competitive immunoassay. In this assay, Standards or unknown urine samples are added to a microplate pre-bounded with E103, and then mixed with anti-E103 monoclonal antibody (mAb). The monoclonal antibody recognizes E103 in Standards or unknown samples under competition in microplate wells. After incubation, the Horse Radish Peroxide (HRP) conjugated secondary antibody is added, followed by an incubation with 3,3',5,5'-tetramethylbenzidine (TMB) substrate. After supplementation of stop solution, the relative reactivity is determined by absorbance measurement at 450 nanometer (nm) wavelength and plotted by comparison with a predetermined E103 standard curve.


**WARNINGS AND PRECAUTIONS**

1. Please note that this is a single-use device. Do not reuse the microplate wells.
2. Do not use reagents after expiration date.
3. Improper temperature exposure during all storage and assay procedures can adversely affect results.
4. Wear protective gloves during all assay procedures.
5. Being light sensitive and skin irritative, TMB Substrate should avoid direct light exposure and skin contact during all storage and assay procedures. In addition, TMB Substrate is flammable, please store it between 2 ~ 8°C and keep it away from any fire source.
6. Avoid skin contact with Stop Solution which contains 0.5 N sulfuric acid. It may cause skin irritation

and burns.








7. Disposal of any discarded materials should be in accordance with local requirements and existing regulations for good laboratory practice.
8. Reliable and reproducible results will be obtained when the assay procedure is practiced with a complete understanding of the package insert instructions and with adherence to good laboratory practice. The test results will be closely related to the operative skills of the end users.
9. The E103 value measured by this product should be corrected before the test result could be evaluated clinically.
10. If experience any serious incident that has occurred in relation to the device, please report to BPM or EC REP as well as the competent authority in your country as soon as possible.

## SYMBOLS GLOSSARY

	Consult Instructions for Use		Do Not Reuse		Use By Date
	Temperature limitation		Catalogue Number		Batch Code
	Date of manufacture		In vitro diagnostic medical device		
	Manufacturer		Authorized representative in the European Community		

## CONTENTS

Sufficient for 96 determinations

COMPONENT	QUANTITY	SYMBOL	Package
Coated Microplate with E103, READY TO USE Handle with care. See Warnings and Precautions.	96 wells: 12 x 8-well strips		Box (-20°C)
Standard, LYOPHILIZED E103 Powder	2 vials		
1000X mAb anti-E103 Monoclonal antibody anti-E103	1 vial, 10 µL		
4000X HRP Conjugate HRP conjugated anti-mouse IgG	1 vial, 10 µL		
Control-H, LYOPHILIZED E103 Powder	2 vials		
Control-L, LYOPHILIZED E103 Powder	2 vials		
Diluent Phosphate buffer with BSA, Tween-20 and ProClin-300	1 bottle, 50 mL		Box (4°C)
10X Wash Buffer Phosphate buffer with Tween-20	1 bottle, 50 mL		
TMB Substrate 0.1% 3,3',5,5'-tetramethylbenzidine with Methanol and Acetone Light sensitive and flammable. See Warnings and Precautions.	1 bottle, 20 mL		
Stop Solution 0.5 N sulfuric acid Handle with care. See Warnings and Precautions.	1 bottle, 20 mL		

## STORAGE AND STABILITY

The shelf life of the DNlite uPTM-FetA ELISA kit is **12 months**, and the components in the reagents should be stored according to the recommendations in the table below.

COMPONENT	Open-vial	
	Expiry date	Temperature (°C)
Standard	Reconstituted Standards (5 µg/mL) are stable for 7 days at -10 ~ -30°C and repeated freeze-thaw cycles should be avoided.	
Control-H & Control-L	Reconstituted Control-H & Control-L are stable for 14 days at -10~-30°C and repeated freeze-thaw cycles should be avoided.	
Coated Microplate	4 weeks (Need to be within the expiry date of each component)	-10 ~ -30°C (Need to be placed in an aluminum foil bag containing desiccant)
1000X mAb anti-E103		-10 ~ -30°C
4000X HRP Conjugate		
Diluent		2 ~ 8°C
10X Wash Buffer		
TMB Substrate		
Stop Solution		

### TRANSPORTATION CONDITION

It is recommended to transport the DNlite uPTM-FetA ELISA Kit at 0°C ~ 9°C for no more than 5 days.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Urinary creatinine assay with regulatory approval or clearance
- Precision single micropipettes or multi-channel pipettes
- Microcentrifuge tubes and disposable tips
- Vortex mixer and microcentrifuge
- Microplate reader capable of endpoint measurement at 450±10 nm
- Plastic container for the preparation of reagents
- Graduated cylinder 500 mL
- Orbital shaker
- Adhesive plate seals
- ddH<sub>2</sub>O

### SPECIMEN COLLECTION AND HANDLING


1. Morning urine samples must be collected in clean and dry containers. Avoid cross contamination.
2. No additives or preservatives are required to preserve the integrity of urine samples.
3. Urine samples can be temporarily stored at 20 ~ 25°C for 3 hours or 2 ~ 8°C for 72 hours before measurement. If sample isn't measured within the 24 to 72-hour period mentioned above, urine samples shall be stored at -10 ~ -30°C for up to 2 weeks.
4. Avoid repeated freezing and thawing of the urine samples.
5. Before performing the assay, please bring the sample temperature back to room temperature.
6. Centrifuge urine samples for 5 minutes at 1,000±20 g. Take supernatant and assay immediately.
7. Use Diluent for sample dilution if necessary.

### REAGENTS PREPARATION

#### A. Prepare Standard Solution

Reconstitute the E103 Standard with 0.2 mL distilled or deionized water, sit for 10 minutes at room temperature until completely dissolved and mix gently. The reconstituted E103 Standard is now at a concentration of 5 µg/mL. Dilute the 5 µg/mL E103 Standard with Diluent by 5-fold, then mix gently

to get a 1 µg/mL E103 Standard. Dilute the 1 µg/mL E103 Standard with Diluent by 2-fold, then mix gently to get a 500 ng/mL E103 Standard. Procedures for the serial dilution to generate Standards (Standard #1 ~ #7) for establishing E103 standard curve are shown in the following table and all Standards must be prepared and mixed well immediately prior to use. Please note that Standard #8 only constitutes 0.2 ml of Diluent.

Content in tube	Conc. (ng/mL)	Standards										
		-	-	-	#1	#2	#3	#4	#5	#6	#7	#8
		5000	1000	500	250	125	62.5	31.25	15.625	7.813	3.906	0
ddH <sub>2</sub> O (mL)		0.2	X	X	X	X	X	X	X	X	X	X
Diluent (mL)		X	0.4	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2

**2 x serial dilution**

Note: (a) Reconstituted Standards (5 µg/mL) are stable for 7 days at -10 ~ -30°C  
(b) Standards #1-#8 should be used after each preparation. Unused portion should be discarded.

### B. Control-H and Control-L

Reconstitute the E103 Controls with 0.2 mL distilled or deionized water, sit for 10 minutes at room temperature until completely dissolved and mix gently.

### C. 1X Wash Buffer

Recover 10X Wash Buffer to room temperature prior to use until all the salt crystals are dissolved. Calculate the required amount of 1X Wash Buffer for each assay. For each microplate, mix 50 mL 10X Wash Buffer with 450 mL distilled or deionized water. Mix uniformly but gently.

### D. 1X mAb anti-E103

Calculate the required amount of 1X mAb anti- E103 for each assay, and mix 1000X mAb anti-E103 with Diluent according to the amount required. For each microplate, mix 8 µL 1000X mAb anti-E103 with 8 mL Diluent. Mix uniformly but gently.

### E. 1X HRP Conjugate

Calculate the required amount of 1X HRP Conjugate for each assay, and mix 4000X HRP Conjugate with Diluent according to the amount required. For each microplate, mix 3 µL 4000X HRP Conjugate with 12 mL Diluent. Mix uniformly but gently.

## ASSAY PROCEDURE

Prepare enough microplate modules for all Standards, Controls, and urine samples and secure them in a holder. **Recover all reagents to room temperature prior to use.**

1. Respectively, add 50 µL of Standard solution #1-#8 (in duplicate), Control-H, Control-L, and urine samples into each well of E103 Coated Microplate.

2. Add 50  $\mu\text{L}$  of 1X mAb anti-E103 into each well with either Standard solution (#1-#8) or urine samples already added.
3. Incubate for 2 hours at 25°C and 200 rpm on an orbital shaker. Add adhesive plate seals to the microplate during incubation to prevent substance from unintentionally falling into the wells.
4. After the incubation, discard the contents in the wells.
5. Wash each well with 300-400  $\mu\text{L}$  1X Wash Buffer. Discard the contents and sharply strike the wells on absorbent paper to remove residual liquid. Wash for a total of 4 times.
6. Add 100  $\mu\text{L}$  1X HRP Conjugate into each well. Incubate the microplate for 30 minutes at 25°C and 200 rpm on an orbital shaker. During incubation add adhesive plate seals as in Step 3.
7. After the incubation, discard the contents in the wells and wash the wells as described in Step 4 and Step 5.
8. Add 100  $\mu\text{L}$  TMB Substrate into each well. Incubate in the dark for 30 minutes at room temperature.
9. Add 100  $\mu\text{L}$  Stop Solution into each well. Mix by a brief shaking until the mixture become homogeneous.
10. Determine the absorbance at  $450\pm 10$  nm within 30 minutes and calculate the results.


## INDICATION OF INSTABILITY OR DETERIORATION OF REAGENTS

When the measurement of the Control-H or Control-L is out of the expected reference range, it may indicate deterioration of the reagents. Associated test results may be invalid and samples need to be retested. Assay recalibration may be necessary.

## DATA ANALYSIS

### A. Calculation

1. The standard curve of this assay fits to a sigmoidal 4- or 5- parameter logistic equation. The results of curve-fitting can be obtained with any computer program having a 4- or 5- parameter logistic function.
2. Graph the standard curve by plotting the OD.450 values on the Y-axis with linear scale against the standard concentrations on the X-axis with logarithm scale. (Use the mean OD.450 values of the duplicate for each concentration of the standard solution.)
3. The concentration of E103 (ng/mL) in patient urine samples can be obtained by inputting the OD.450 value and inversing fitted standard curve.

**B. Quality criteria of assay**

The quality control (Control-H and Control-L) should be included in each experiment for quality assurance. The suggested reference range of quality control is listed below. (Each testing lab shall establish its own reference range.)

	Control-H	Control-L
E103 value (ng/mL)	82.78 ~ 144.35	11.55 ~ 30.87

**C. Clinical application**

All concentration of E103 need to be urinary creatinine-corrected for clinical use. The calculation is as follows:

$$E103/UCr = \frac{E103 \text{ (ng/mL)}}{\text{urine creatinine (mg/mL)}}$$

**D. The risk of kidney function decline progression**

1. The risk classification of kidney function decline progression in **type 2 diabetes patients** shall be based on E103/UCr. Please follow the procedure below:
  - I. If E103(ng/mL) < 5.428 ng/mL (LoQ), then E103 shall be reported as “<5.428 ng/mL”, and classify as “low risk”.
  - II. If E103(ng/mL) ≥ 250 ng/mL (upper limit of measuring range), then E103 shall be reported as “≥ 250 ng/mL” and classify as “high risk”.
  - III. If E103(ng/mL) is between 5.428 ~ 250 ng/mL, then calculate E103/UCr ◦
    - i. If E103/UCr < 7.53 ng/mg, then classify as “low risk”.
    - ii. If E103/UCr ≥ 7.53 ng/mg, then classify as “high risk”.
  
2. The risk classification of **graft failure in kidney transplantation recipients** shall be based on E103/UCr. Please follow the procedure below:
  - I. If E103(ng/mL) < 5.428 ng/mL (LoQ), then E103 shall be reported as “< 5.428 ng/mL”, and classify as “low risk”.
  - II. If E103(ng/mL) ≥ 250 ng/mL (upper limit of measuring range), then E103 shall be reported as “≥ 250 ng/mL” and classify as “high risk”.
  - III. If E103(ng/mL) is between 5.428 ~ 250 ng/mL, then calculate E103/UCr ◦
    - i. If E103/UCr < 29.1 ng/mg, then classify as “low risk”.
    - ii. If E103/UCr ≥ 29.1 ng/mg, then classify as “high risk”.

**ANALYTICAL PERFORMANCE**

**1. Detection capability**

Limit of Blank (LoB) = 0.216 ng/mL;  
 Limit of Detection (LoD) = 1.901 ng/mL;  
 Limit of Quantitation (LoQ) = 5.428 ng/mL

**2. Linearity**

4.90 ~ 255.39 ng/mL

**3. Measuring range**

5.428 ~ 250 ng/mL

#### 4. Precision

##### i. Repeatability

Samples were tested in triplicate. Reagents from three different lots were used each day for 20 days. Variation of within-run, run-to-run, and day-to-day are listed below.

Sample	Mean conc. of E103 (ng/mL)	Within-run (%CV)	Run-to-run (%CV)	Day-to-day (%CV)	Total CV
Sample A	11.32	8.12	8.45	6.69	13.49
Sample B	32.09	7.57	11.63	0	13.88
Sample C	40.46	4.06	9.1	0	9.97
Sample D	59.39	5.22	3.42	4.45	7.67
Sample E	130.38	5.2	4.31	4.6	8.17
Sample F	172.92	5.51	9.17	0	10.69

##### ii. Reproducibility

Samples were tested in triplicate with two runs each day. Tests were performed in three independent laboratories for 5 days. Variation of within-run, run-to-run, day-to-day, and site-to-site are listed below.

Sample	Mean conc. of E103 (ng/mL)	Within-run (%CV)	Run-to-run (%CV)	Day-to-day (%CV)	Site-to-Site (%CV)	Total CV
Sample a	8.72	12.38	9.13	16.95	10	24.98
Sample b	26.83	8.98	6.64	13.74	13.91	22.51
Sample c	39.71	8.6	5.56	12.55	15.17	22.2
Sample d	67.26	6.34	5.04	10.62	11.27	17.48
Sample e	142.57	5.66	5.11	8.33	11.05	15.8
Sample f	167.38	6.13	4.91	9.01	10.81	16.12

#### 5. Interference

Potential interference	The highest non-interference concentration (mg/L)
Acetaminophen	600
Ascorbic Acid	90
Albumin	1500
Direct bilirubin	1800
Creatinine	3750
Enalapril	187.5
Glucose	60000
Glibenclamide	6
Hemoglobin	375
Ibuprofen	750
Losartan	300
Metformin HCL	20000
Metronidazole	120
Salicylic Acid	37.5
Simvastatin	45
Sodium Citrate	150
Sodium oxalate	750
Trichloromethiazide	75
Urea	6250
Uric Acid	375
Urobilinogen	37.5

## 6. Metrological Traceability

The E103 Standard value is assigned to be traceable to arbitrarily defined reference materials (ARM) according to international standard ISO17511. The ARM is an analyte manufactured by BPM. The relative bias of all batches E103 standard are less than 15 %. The E103 Control-H and Control-L values are assigned to be traceable to E103 Standard.

## CLINICAL PERFORMANCE AND APPLICATION

These clinical performance studies evaluated the use of E103/UCr (ng/mg) as a marker for (1) kidney disease progression in type 2 diabetes patients and (2) graft failure in kidney transplantation recipients (KTRs).

For (1), the clinical cut-off for E103/UCr was determined in a primary cohort of 294 Taiwanese patients (mean age 61, 55% male) as 7.53 ng/mg. Results showed that patients with E103/UCr above 7.53 ng/mg had a higher risk of kidney function decline compared to those with low values (hazard ratio 8.94, p-value<0.001). Results were validated in a separate cohort of 376 patients from the Netherlands (mean age 64, 66% male) where similar findings were observed. High-risk patients had a significantly higher risk of worsening renal function (hazard ratio 1.87, p-value <0.001).


For (2), the E103/UCr clinical cut-off 29.1 ng/mg was determined in a primary set of 267 KTRs (mean age 53.9, 57% male). Results showed that KTRs with E103/UCr above 29.1 ng/mg had a higher risk of kidney graft failure than those with low values (hazard ratio 4.98, p-value<0.001). In the validation set of 267 KTRs (mean age 52.32, 58% male), KTRs with E103/UCr above 29.1 ng/mg had a significantly higher risk of kidney graft failure than those with low E103/UCr values (hazard ratio 4.51, p-value <0.001).

## REFERENCES

1. Chuang, GT et al. (2023) Urinary Fetuin-A Fragments Predict Progressive eGFR Decline in Two Independent Type 2 Diabetes Cohorts of Different Ethnicities. Am J Nephrol, DOI: 10.1159/000534514
2. Alkaff, FF et al. (2023) Urinary Post-Translationally Modified Fetuin-A Protein Is Associated with Increased Risk of Graft Failure in Kidney Transplant Recipients. Am J Nephrol, DOI: 10.1159/000534829

**Bio Preventive Medicine Corp.**  
6F-1, 6F-2, 6F-3, 6F-5, 6F-6,  
No.18 Taiyuan St., Zhubei City,  
Hsinchu County, 30265 Taiwan

**EMERGO EUROPE**  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands

### AUSTRALIAN SPONSOR:

**Emergo Australia**  
Level 20 Tower II, Darling Park  
201 Sussex Street  
Sydney, NSW 2000, Australia

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